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Instructions for Use

Repair and customization:

PMMA for brain can be repaired or customized with the following products:

- PMMA for brain Dentin/Enamel: Tooth-coloured PMMA plastic
- PMMA for brain Dentin Liquid
- PMMA for brain Bonder: Adhesive

Contraindications/interactions:

In patient hypersensitivity to one of the ingredients, these products may not be used at all or only under the particular scrutiny of the dentist or physician in charge. The composition is available on request if needed. Known cross-reactions or interactions of this medical product with other materials already present in the oral environment must be taken into consideration by the dentist in charge when using this product.

Please note the safety instructions on the back of these Instructions for Use!

The product properties are based on adherence to and compliance with these Instructions for Use.

Processing notes:

PMMA for brain Dentin/Enamel

Dosage:

Mixing ratio PMMA for brain Dentin/Enamel: Liquid, 1 g : 0.7 ml

Mixing time 30 s*

Swelling time 15 s to max. 45 s following the mixing time*

Flow phase 240 s

Polymerization 15 min in a pressure pot at 45 °C and 2–2.5 bar pressure

*All values assume a room temperature of 21 °C.

Wherever possible, use class IV dental stone for cast fabrication. Always wet the cast thoroughly and isolate the plaster from the using a commercially available alginate.

Bonding:

- Degrease and abrade the contact surfaces before processing.
- Avoid any contamination by wax or isolating compounds.

PMMA for brain Bonder:

- Remove residual wax, grease, grinding dust and other contaminants before applying. After steam-cleaning ensure that all parts, including undercut, are thoroughly dry. When drying with compressed air, make sure that no contamination occurs (e.g. by oil residue from the compressed-air system).
- Roughening the surfaces to be joined. For larger surfaces, it may be advisable to provide small retention grooves.

- Wet the surface to be bonded with PMMA for brain Bonder and apply enough bonder for the contact surface to stay moist for 5 min. If the contact surface has dried of all PMMA for brain Bonder during this time, re-moisten with PMMA for brain Bonder. The exposure time must not be less than 5 min. PMMA for brain Dentin/Enamel must be applied no later than 10 minutes after conditioning the surface.

- Avoid any contact with plaster or with materials containing silicone.

Finishing and polishing:

- To avoid inaccuracies of fit after polymerization, avoid excessive heat during finishing and polishing.

Special notes:

- Close containers tightly after use.
- Do not use the material after the stated expiration date.
- Store in a dry and cool place (up to 25 °C). Avoid direct sunlight.
- Indicate the batch (LOT) number for every transaction that requires an identification of the material.
- Notify the dentist or physician in charge of all factors described above if you use this medical product for a custom construction. When working with these materials, make sure to comply with the product-specific Material Safety Data Sheets (MSDS).

For detailed information, please download the relevant Safety Data Sheet from our website download area at www.degudent.com

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